

RESIN 665

NEW

Cold setting resin for sample recovery.

Mounting is a step in metallographic sample preparation which enables samples to be ground and polished more ergonomically. Embedding the sample in a resin is usually irreversible and recovering the sample (if necessary) after embedding proves almost impossible.

Description

Resin 665 is a cold mounting resin that has the unique ability to dissolve in acetone. This ability is transformed into a sustainable functionality to recover metallographic specimens after mounting. After metallographic processing of costly materials such as noble metals, it is often important from an economic and technical perspective to recover the samples from the mounting resin.

Characteristics

Resin 665 is an easy to use acrylic bi-component resin (powder and liquid). During curing, the mixture can reach up to 105°C (embedding mould 30 mm). The hardness of the cured resin is about 80 Shore D. Like LAM PLAN's other acrylic cold setting resins (Resin 605 and Resin 609), Resin 665 can be used to embed all types of material for metallographic analysis.

Recovering the sample


There are two ways to recover the specimens from the mount:

1. Using Acetone:

The sample can be placed into a closed receptacle containing acetone at ambient temperature (do not heat the acetone). Make sure that the sample is totally submerged in the liquid. (approximately 150 ml of acetone is recommended to dissolve a 30 mm mount). Depending on the size of the mount, it can take from a few hours to three days for the embedding resin to dissolve.

2. Using heat:

This method can only be used if the specimens are resistant to temperatures above 150°C. The embedded sample can be placed in an oven preheated to 150°C. In a few hours depending on the size of the mount, the resin softens allowing the sample to be ejected mechanically.



Technical Data	RESIN 665
Packaging	Kit powder 1000 g + catalyst 500 ml Code 06 00665 00
	Powder 1000 g - Code 06 01665 00
	Catalyst 500 ml - Code 06 02665 00
Dosage	2 parts powder to 1 part liquid
Polymerisation time	10 to 15 min.
Polymerisation temperature	Maximum 105°C
Hardness	80 Shore D
Dissolution	Acetone
Demoulding temperature	150°C



Safety data sheets