crosscopy, the support ring is manufactured with a ceramic surface parallelism of ± 0.02 mm.

The ceramic material is extremely hard and cannot be machined using silicon carbide (SiC) paper. If correction of dimensional accuracy becomes necessary, the ceramic surface can be hand-lapped using a diamond wheel or a lapping plate.

The manual sample holder is not intended for use on diamond wheels, lapping wheels, grinding machines, or similar equipment.


The required thin section/sample should be prepared manually on SiC paper under water cooling. The specimen thickness can be adjusted by means of the ceramic ring. One scale division corresponds to 0.02 mm.

Scope of Supply - For Metallography

- Hand specimen holder, type M.
- Inserts for specimen diameters $\varnothing$ 25 mm and $\varnothing$ 30 mm.

For Lithopreparation

- Hand specimen holder, type H.
 - Insert for specimen diameter $\varnothing$ 40 mm.
- Special versions are available upon request.

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MANUAL SAMPLES HOLDER

METALLOGRAPHY - LITHOPREPARATION

1. Loosen the set ring and turn the setting ring forward until the scale reading is zero and the sample locking screw, operated with a 5 mm Allen key, appears in the ring bore.
2. Lock this position by tightening the set ring.
3. Place the sample holder, with the ceramic ring, on a flat surface. Position the sample in the sample holder so that the surface to be prepared rests on a flat surface.
4. Press the sample holder and sample on to the flat surface, then secure the sample by tightening the sample locking screw.
5. Loosen the set ring and turn the setting ring back to the required grinding depth. Lock the setting by tightening the set ring.

If the ceramic ring becomes uneven due to frequent use, it can be restored to a flat condition by flat grinding.

Inserts for the Manual Sample Holder

1. Unscrew the setting ring from the manual sample holder and remove the sample locking screw. Turn the set ring back until the locking screw becomes visible.
2. Push the insert into the sample holder and rotate it until the locking screw of the insert sleeve becomes visible through the threaded opening of the holder. Secure the insert sleeve with the set screw.
3. Refit the setting ring to the manual sample holder.
4. Operate the manual sample holder as described above.

When installing inserts, push them in gently and keep them aligned. Do not tilt them. The insert is manufactured to an H7 fit and should slide into place without force.

Cleaning

After use, the manual sample holder should be cleaned and the threads regreased if necessary. Dirt or abrasive residue in the threads may cause the setting ring to stick, making the manual sample holder unusable.

Note

Repairs resulting from damage to the ceramic ring cannot be carried out free of charge. Damage to the ceramic ring is generally caused by improper use.

As the name indicates, this is a manual sample holder and is intended to be operated manually. It is designed for the preparation of thin sections mounted on glass slides, for example in histology for the microscopic examination of bone tissue.

To ensure the required parallelism of thin sections for light mi-